

# STERILSEAL

automatic thermosealer



## STERILSEAL

	M.U.	
Packet tray		optional
Overall dimensions W x H x D	mm	380 x 100 x 240
Electrical voltage and frequency	V - Hz	230 - 50/60
Absorbed maximum power	W	800
Seal width	mm	11
Seal duration s 3	m/min	10
Weight (mass)	kg	7

Technical data are subject to change without notice



STERILSEAL is the automatic sealer with adjustable temperature (temperature range: 50°C - 180°C). STERILSEAL is the ideal solution for carrying out many weldings per day (time can be saved). The ergonomic entirely stainless steel hull design guarantees a simple use. It carries out a fast sealing procedure and can be equipped with an optional packing support.